# INHEALTH<sup>®</sup>

## TECHNOLOGIES

### Blom-Singer<sup>®</sup> Plug Insert

USER Instructions For Use

### Blom-Singer<sup>®</sup>

voice restoration systems

R1

37582-01G

37582-01G Effective October 2017 / Gültig ab Oktober 2017 / Efectivo a partir de octubre de 2017 / Date d'entrée en vigueur : octobre 2017 / Valido da ottobre 2017 / Met ingang van oktober 2017 / Vigente em outubro de 2017 / V platnosti od října 2017 / Wchodzi w życie w październiku 2017 r.

**Blom-Singer** and **InHealth Technologies** are registered trademarks of **Freudenberg Medical, LLC**

# TABLE OF CONTENTS

Inhaltsverzeichnis / Índice / Sommaire / Sommario / Inhoud / Obsah / Spis treści

ENG	<b>Blom-Singer Plug Insert</b>	<b>4</b>
DEU	<b>Blom-Singer Stopfeneinsatz</b>	<b>12</b>
SPA	<b>Obturador Blom-Singer</b>	<b>22</b>
FRA	<b>Insert d'occlusion Blom-Singer</b>	<b>31</b>
ITA	<b>Tappo Blom-Singer</b>	<b>40</b>
NLD	<b>Blom-Singer-Insteekstop</b>	<b>49</b>
POR	<b>Suplemento com Plugues Blom-Singer</b>	<b>58</b>
CES	<b>Zátková Vložka Blom-Singer</b>	<b>67</b>
POL	<b>Wtyk Blom-Singer</b>	<b>75</b>

## **Blom-Singer®**

voice restoration systems

### **USER Instructions For Use**

Gebrauchsanweisung für BENUTZER  
Instrucciones para el USUARIO  
Instructions d'utilisation destinées à L'UTILISATEUR  
Istruzioni per l'uso destinate ALL'UTENTE  
Gebruiksaanwijzingen voor de GEBRUIKER  
Instruções de uso para o USUÁRIO  
Návod k použití pro UŽIVATELE  
Instrukcja obsługi dla UŻYTKOWNIKA

### **BLOM-SINGER® PLUG INSERT**

#### **Introduction**

*Please refer to the diagrams located at the back of this instruction manual. A glossary of terms appears at the end of this chapter.*

**These products are intended for single patient use only. Not for use between patients.** See Warnings, Precautions and Complications.

#### **INDICATIONS**

The Blom-Singer Plug Insert may be used on a temporary basis in conjunction with a Blom-Singer indwelling voice prosthesis to prevent leakage of fluids through the voice prosthesis into the trachea during swallowing.

#### **CONTRAINDICATIONS**

The Plug Insert is a medical product and should be used only by individuals with experience and training in its use and care.

The Plug Insert is not intended to permanently replace a properly functioning voice prosthesis. It is not intended to be placed directly in a surgically created fistula between the esophagus and trachea. It is not intended to be used in combination with any other voice prosthesis or device other than the Blom-Singer ADVANTAGE or CLASSIC Indwelling Voice Prostheses.

**Warning: Use of all tracheostoma accessories is contraindicated (must not be used) when the Plug Insert is placed within the voice prosthesis. Use of tracheostoma accessories in conjunction with the Plug Insert may result in accidental displacement of the device or voice prosthesis, which could result in inhalation of the device or voice prosthesis.**

**Do not use the Plug Insert with a 4 mm size voice prosthesis. The Plug Insert will not fully insert into a 4 mm size voice prosthesis. This could result in pushing the voice prosthesis into the esophagus.**

## **HOW SUPPLIED**

The Blom-Singer Plug Insert is supplied non-sterile as an accessory device for use with a Blom-Singer indwelling voice prosthesis. An inserter tool is included with each device.

## **PRODUCT DESCRIPTION**

The Blom-Singer Plug Insert is made of silicone, a biocompatible material used in many medical products. It is available in two models (16 Fr. and 20 Fr.) for use with the Blom-Singer Classic and Advantage Indwelling Voice Prostheses.

The Plug Insert (diagram 1) is designed to prevent leakage through the voice prosthesis until such time that the voice prosthesis and TEP is evaluated by a clinician; it does not allow voicing. The Plug Insert consists of a plug tip (diagram 1a) and neck strap with safety hole (diagram 1b). The Plug Insert may be inserted by the user.

## **WARNINGS AND PRECAUTIONS**

The qualified, trained medical professional is to instruct patients on the use and care of this device and to provide patients with the instructions for use supplied with the product.

Should the device show signs of malfunction, contact your physician/clinician as soon as possible to have the device performance evaluated.

If there are tears, cracks, or structural damage to the device, discontinue use.

Changes in your anatomy or medical status may lead to improper fitting and/or function of the device. It is recommended that your voice prosthesis and TEP be evaluated by your clinician at regular intervals.

This device is designed for temporary use only (maximum of 30 days) and is not intended to replace the indwelling voice prosthesis. After placement of the device, contact your physician/clinician as soon as possible.

Do not detach or in any way modify the neck strap on the Plug Insert.

Do not use solvents or petroleum-based products for cleaning or lubricating the device. These materials may damage silicone or

---

cause the device to not work properly.

Always clean the device after it has been removed per instructions provided in the cleaning section below.

Only use a lint free cloth or tissue to dry the device. Use of non-lint free materials may leave debris which can be aspirated into the airway.

Never force the Plug Insert into the voice prosthesis. If the device does not insert easily on the first attempt, do not continue to try to insert it. Should you experience repeated difficulty in managing the Plug Insert, contact your physician/clinician as soon as possible.

It is very important to not insert the Plug beyond this point of palpable resistance. This could result in pushing the voice prosthesis into the esophagus.

### **Voice Prosthesis Dislodgement**

Care should be exercised when inserting or removing the Plug Insert, or inserter tool to avoid accidental displacement of the voice prosthesis which could result in aspiration (inhalation) of the prosthesis. In the unlikely event that this should occur, you should attempt to cough the voice prosthesis out of the windpipe. Immediately seek medical attention if unsuccessful.

The neck strap should not, under any circumstances, be cut or removed from the Plug Insert. The neck strap is intended to help prevent accidental aspiration (inhalation) of the device should the device become dislodged.

If the voice prosthesis is accidentally dislodged from the TE puncture, immediately place a Blom-Singer Puncture Dilator, or a suitable device as recommended by your doctor, of the appropriate diameter, in the puncture to keep it from closing and leaking fluids. Contact your doctor immediately. The voice prosthesis should be reinserted within 24 hours.

Never insert foreign objects into the voice prosthesis. Inserting objects other than the supplied device may cause dislodgment of the voice prosthesis or its components and may cause these objects to be swallowed or inhaled.

Never insert any cleaning device while the Plug Insert is in place.

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

## Complications

Although rare, the following complications have been identified to occur with Blom-Singer Plug Inserts. They include:

- stoma or device contamination or bacterial infection, which may require removal of the prosthesis and or treatment with appropriate antibiotics;
- accidental aspiration of the voice prosthesis into the airway, which may require removal by a physician;
- aspiration of liquids through the prosthesis into the airway, which may cause the patient to cough;
- occasional movement out of position by the voice prosthesis, requiring replacement after enlargement of the TE puncture and additional supervision of the stoma care regimen;
- dislodgment of the voice prosthesis and subsequent closure of the TE puncture;
- tearing or other damage to the voice prosthesis from improper use;
- accidental swallowing of the voice prosthesis into the esophagus and/or gastrointestinal (GI) tract.

## INSTRUCTIONS FOR USE

Your clinician will instruct you on the routine use and care of this Blom-Singer device.

Changes in your anatomy or medical status may lead to improper fitting and/or function of the device. It is recommended that your voice prosthesis and TEP be evaluated by your clinician at regular intervals.

If there are tears, cracks, or structural damage to the device, discontinue use.

The indwelling accessory device **should always be** cleaned to prevent contaminants from entering into the tracheostoma, puncture, or esophagus, see below for cleaning instructions.

The user's hands should always be washed with soap and water prior to cleaning the Plug Insert.

---

*The following procedural instructions are provided by Eric D. Blom, Ph.D. Please refer to the diagrams located at the back of this instruction manual.*

### **Insertion of the Plug Insert**

**Ensure that the size of the Plug Insert correctly matches size of the voice prosthesis.**

Insert the Plug Insert only while using a close-up mirror and a bright light focused directly into the stoma. The user's hands should always be washed with soap and water prior to cleaning the plug insert.

1. Place the tip of the inserter tool into the open end of the Plug Insert (diagram 2). Handle the device by the base of the strap to avoid touching the portion that inserts into the voice prosthesis.
2. Attach (diagram 3) the safety hole in the neck strap on the Plug Insert securely and completely down over the safety peg on the inserter tool. This should prevent inadvertent dislodgment of the Plug Insert off the inserter tool during insertion into the voice prosthesis.
3. Hold the end of the inserter tool as shown in diagram 4 and align the top of the Plug Insert partially inside the voice prosthesis with the neck strap oriented upwards (diagram 5).
4. Carefully and gently insert the Plug Insert until it is fully inserted and positioned flush with the voice prosthesis (diagram 6).

Never force the Plug Insert into the voice prosthesis.

**Caution:** If the device does not insert easily on the first attempt, do not continue to try to insert it. Should you experience repeated difficulty in managing the Plug Insert, contact your physician/clinician as soon as possible.

**Caution:** Do not use Plug insert with a 4 mm size voice prosthesis. The Plug Insert will not fully insert into a 4mm size voice prosthesis due to the shorter length of the voice prosthesis shaft. This could result in pushing the voice prosthesis into the esophagus.

5. After confirming that the Plug Insert is securely in place and fully inserted, the inserter tool should be carefully removed.
6. Hold the inserter tool and the neck strap together firmly in one hand. Carefully lift the neck strap off the safety peg on the inserter tool. Hold the neck strap against the neck while slowly and gently withdrawing the inserter tool with a twisting motion (diagram 7).

**Warning:** Caution should be exercised when removing the inserter tool to avoid accidental displacement of the Plug Insert or voice prosthesis which could result in aspiration (inhalation) of the device or voice prosthesis. Should aspiration occur, attempt to cough the device and/or voice prosthesis out of the windpipe. Seek medical attention **immediately** if unsuccessful.

7. Apply a strip of hypo-allergenic, medical adhesive tape over the neck strap to prevent strap movement as shown in diagram 8.

**Warning:** Do not detach or in any way modify the neck strap on the Plug Insert. **The neck strap should not, under any circumstances, be cut or removed from the Plug Insert.** The neck strap is intended to help prevent accidental aspiration (inhalation) of the device should the device become dislodged.

### **Removal of the Plug Insert**

Removal of the Plug Insert should only be done while using a close-up mirror and a bright light focused directly into the stoma so that the open end of the prosthesis is clearly visible. The users hands should always be washed with soap and water prior to cleaning the Plug Insert.

1. Stabilize the voice prosthesis by holding it carefully and firmly with tweezers to avoid displacement (diagram 9).
2. Grasp the neck strap of the Plug Insert and gently pull it out from the voice prosthesis.

## **CLEANING AND CARE INSTRUCTIONS**

Your clinician will instruct you on the routine use and care of this Blom-Singer indwelling accessory device.

---

After use, the Plug Insert should be inspected and any visible debris should be removed. Remove visible debris by holding the Plug Insert under a running stream of water.

The user's hands should always be washed with soap and water prior to cleaning the Plug Insert.

Only use a lint free cloth or tissue to dry the device. Use of non-lint free materials may leave debris which can be aspirated into the airway.

Once the Plug Insert has been cleaned and dried completely, it should be stored in a clean, resealable plastic bag. Please refer to the instructions for use supplied with the indwelling voice prosthesis for complete instructions on proper care and maintenance of the voice prosthesis.

Dispose of the Plug Insert (non-biodegradable material) in accordance with all federal, state, and local regulations. Responsibility for proper waste disposal is with the owner of the waste.

## **ORDERING INFORMATION**

### **USA**

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337, Monday–Friday, 9:30am–7:00pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: [order@inhealth.com](mailto:order@inhealth.com) ORDER ON-LINE: [www.inhealth.com](http://www.inhealth.com) POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

### **International**

Please contact our customer service department for a distributor referral.

### **Consumer Affairs**

If you have any questions or dissatisfaction with a product please contact our customer service department via telephone, fax, post, or Email: [productcomplaints@inhealth.com](mailto:productcomplaints@inhealth.com)

## GLOSSARY

### **Airway**

The air passage by which air reaches the lungs.

### **Clinician**

A health professional, such as a physician, licensed speech-language pathologist, or nurse who is involved in clinical practice.

### **Contraindications**

A reason that makes it inadvisable to prescribe a particular procedure or treatment.

### **Esophagus**

The muscular membranous tube for the passage of food from the pharynx to the stomach; the “food pipe.”

### **Fistula; TE-fistula**

An abnormal duct or passage that connects an abscess or cavity to the body surface; a small opening between the trachea and the esophagus.

### **Indications**

A reason to prescribe a drug or perform a procedure.

### **Tracheoesophageal puncture (TEP)**

An opening made by a surgeon between the trachea and the esophagus. The opening is made at the time of initial surgery or later so that a voice prosthesis can be placed into this passage.

### **Trachea**

The airway through which respiratory air passes to and from the lungs (windpipe).

### **Tracheoesophageal (TE)**

Pertaining to the trachea and esophagus.

### **Tracheostoma**

An opening through the neck and into the trachea. It is through this opening that the laryngectomee breathes.

### **Voice prosthesis**

A one-way speech valve that is inserted into the TEP to allow air to pass from the trachea through the prosthesis and into the esophagus to produce sound.

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