

# Nasal Septal Perforation Prosthesis

**Blom-Singer®**

The voice of experience since 1978

*The new solution for  
non-surgical closure*

## 2-piece magnetically coupled prosthesis

A study<sup>1</sup> of patients using a 2-piece magnetically coupled prosthesis showed:

- *No particular discomfort on insertion, removal, or replacement*
- *Easier debris removal for ENTs*

## Choice of shape, size, and thickness for a comfortable fit

Two medical-grade, silicone-covered neodymium magnets couple securely.

Multiple options – round or oval, 6 sizes each – help ensure an accurate fit:

- *Less pistoning*
- *Less side-to-side movement*

## Reduce or eliminate symptoms:

- *Bleeding, crusting*
- *Breathing difficulty, nasal pressure, whistling*
- *Rhinorrhea, epistaxis, infection*



# Nasal Septal Perforation Prosthesis


The new solution for non-surgical closure

**Blom-Singer®**


The voice of experience since 1978

## OVAL



### 11 x 13 mm



### 13 x 17 mm



### 15 x 21 mm




Flange width is 5mm for all prostheses.


**Diagrams actual size**

## ROUND


### 5 mm



### 7 mm



### 10 mm



**INHEALTH**  
TECHNOLOGIES®

**severn**   
SEVERN HEALTHCARE TECHNOLOGIES

1. Teschner M, et al. *Eur Arch Otorhinolaryngol.* 2012 Mar;269(3):861-5.

The Blom-Singer Nasal Septal Perforation Prosthesis was designed by Eric D. Blom, Ph.D., co-inventor of the breakthrough in voice restoration: the Blom-Singer Voice Prosthesis.

©2017 InHealth Technologies — Manufactured by Freudenberg Medical, LLC (161011)

Distributed in the UK and Ireland by  
Severn Healthcare Technologies  
[www.severnhealthcare.com](http://www.severnhealthcare.com)

Email: [service@severnhealthcare.com](mailto:service@severnhealthcare.com) Tel:  
01635 887640